

```
In [ ]: from PIL import Image
        from IPython.display import display
```

# Single Imputer Report

## Overview

Completed several operations working on single daphnia tracking:

- Created a class `NPZer` which unzips and unpacks `.npz` files, the output from TRex, for imputation and transformation.
- Created a class `TRexDataCleaner` which removes disjoint data points, such as large jumps in position and missing data.
- Created a class `TRexImputer` which provides a framework for applying different imputation strategies to missing daphnia tracking data, currently using only an `avgValue` function.
- Developed `plotDetail` : an easy to use and standardized plotting function for displaying daphnia tracking.
- Established `avgVelocity` a framework for estimating maximum daphnia velocity used in data cleaning and imputation.

---

## Walkthrough of Imputation Process

This process occurs *POST* TRex animal tracking for a single daphnia:

### 1. Unzipping and Preparation

- Use a standardized unzipping function.
- Include functionality to flip axes for plotting if necessary.

### 2. Data Cleaning

- Remove discontinuities in the data.
- Focus on critical clean-up, such as removing data where the tracker has jumped to the edge of the dish.

### 3. Imputation


- Apply an imputation function with modular strategies.
- Ensure the function can accommodate various imputation methods as they are developed.

### 4. Standardized Plotting

- Implement consistent plotting of the daphnia tracking data.

- Ensure plots accurately reflect cleaned and imputed data.

```
In [ ]: from PIL import Image
from IPython.display import display
#display output output\email_update_single_imputer\output_images\side_by_side_imput
display(Image.open('output_images/side_by_side_imputation.png'))
```



## NPZer.py

`NPZer.py` contains the class `NPZer`, which unzips a `.npz` file (output of TRex tracking). This data can be directly converted to a pandas table and subsequently to a `.csv` file or plotted.

### Example Usage:

```
# Import necessary tools
import pandas as pd
import numpy as np
from src.data_manipulation.NPZer import NPZer

NPZer = NPZer()

# Set desired parameters
SOURCE_DIR = 'data/clean_fish_data/fish_data_clean.csv'
INVERT_Y = True
PARAMS = ['time', 'X', 'Y']

# Unzip and turn data into a pandas table
unzippedData = NPZer.pandafy(source_dir=SOURCE_DIR, invertY=INVERT_Y,
params=PARAMS)

# Print data in the form of a pandas table
print('TRex Data:\n', unzippedData)
```

## Output

```
In [ ]: display(Image.open('output_images/NPZer_output.png'))
```

```
TRex Data:
      time      X      Y
0    0.000000  23.536650 -1.792803
1    0.016949  23.517750 -1.792841
2    0.033898  23.517750 -1.792841
3    0.050847  23.517750 -1.792841
4    0.067796  23.517750 -1.792841
...
10817  183.338989  19.579285 -6.965172
10818  183.355927  19.569004 -6.989434
10819  183.372879  19.588287 -7.017863
10820  183.389832  19.577187 -7.083682
10821  183.406784  19.577187 -7.083682

[10822 rows x 3 columns]
```

## TRexDataCleaner.py

TRexDataCleaner.py contains the class TRexDataCleaner, which is used to clean disjoint data in TRex tracking data for later imputation. Disjoint data points are rendered as np.inf.

Example Usage:

```
# Import necessary tools
from src.data_manipulation.TRexDataCleaner import TRexDataCleaner

dataCleaner = TRexDataCleaner()

# Set desired parameters
VMAX = 1

# Set sample of original data
originalData = unzippedData[:25]
```

```

# Print sample of original data
print('Original Data:\n', originalData)

# Clean data
cleanedData, removedData =
dataCleaner.renderDiscontinuities(data=originalData, vmax=25)

# Print cleaned data
print('Cleaned Data:\n', cleanedData)

```

---

## TRexImputer.py

TRexImputer.py holds the class TRexImputer which is used to impute/fill in discontinuous TRex data with a desired impute function. This should be used with data already cleaned with TRexDataCleaner

## Example Usage

```

# Import necessary tools
from src.data_manipulation.TRexImputer import TRexImputer

imputer = TRexImputer()

# Set desired parameters
DATA = cleanedData
FUNCTION = 'avgValue'

# Print original data
print('Original Data:\n', cleanedData)

# Impute data
imputedData = imputer.impute(data=DATA, function=FUNCTION)

# Print imputed data
print('Imputed Data:\n', imputedData)

```

```
In [ ]: display(Image.open('output_images/TRexImputer_output.png'))
```

## Original Data:

	time	X	Y
0	0.000000	23.536650	-1.792803
1	0.016949	inf	inf
2	0.033898	23.517750	-1.792841
3	0.050847	23.517750	-1.792841
4	0.067796	23.517750	-1.792841
5	0.084745	23.529652	-1.792807
6	0.101694	23.525852	-1.792829
7	0.118644	23.525852	-1.792829
8	0.135593	23.525852	-1.792829
9	0.152542	23.529652	-1.792807
10	0.169491	23.525852	-1.792829
11	0.186440	23.529652	-1.792807
12	0.203389	23.529652	-1.792807
13	0.220338	23.529652	-1.792807
14	0.237288	23.529652	-1.792807
15	0.254237	inf	inf
16	0.271186	23.529652	-1.792807
17	0.288135	inf	inf
18	0.305084	23.529652	-1.792807
19	0.322033	23.529652	-1.792807
20	0.338983	inf	inf
21	0.355932	23.525852	-1.792829
22	0.372881	23.529652	-1.792807
23	0.389830	23.536650	-1.792803
24	0.406779	23.540424	-1.792785

Imputing with: avgValue()

## Imputed Data:

	time	X	Y
0	0.000000	23.536650	-1.792803
1	0.016949	23.527200	-1.792822
2	0.033898	23.517750	-1.792841
3	0.050847	23.517750	-1.792841

4	0.067796	23.517750	-1.792841
5	0.084745	23.529652	-1.792807
6	0.101694	23.525852	-1.792829
7	0.118644	23.525852	-1.792829
8	0.135593	23.525852	-1.792829
9	0.152542	23.529652	-1.792807
10	0.169491	23.525852	-1.792829
11	0.186440	23.529652	-1.792807
12	0.203389	23.529652	-1.792807
13	0.220338	23.529652	-1.792807
14	0.237288	23.529652	-1.792807
15	0.254237	23.529652	-1.792807
16	0.271186	23.529652	-1.792807
17	0.288135	23.529652	-1.792807
18	0.305084	23.529652	-1.792807
19	0.322033	23.529652	-1.792807
20	0.338983	23.491658	-1.793025
21	0.355932	23.525852	-1.792829
22	0.372881	23.529652	-1.792807
23	0.389830	23.536650	-1.792803
24	0.406779	23.540424	-1.792785

### avgValue.py

avgValue.py holds the `impute` function which imputes data in between gaps of discontinuous points. The function calculates the velocity between each gap of data, and imputes the respective points according to it.

### Example usage

Usage is shown in above example as `FUNCTION = 'avgValue'`

---

### plotDetail.py

plotDetail.py contains the `plotDetail` function, which:



- Accepts a CSV or NPZ file.
- Requires declarations of X, Y, and time values.
- Plots a graph representing the path of the Daphnia's movements with placeholder labels.

## Example Usage

```
# Import necessary tools
from src.data_manipulation.plotDetail import plotDetail
import pandas as pd
from matplotlib import pyplot as plt
import numpy as np

# Copy the path of the NPZ file
npz_path =
r"/Users/ibrahimrahat/Documents/GitHub/daphnia/data/npz_file/single_7_9_fish1

# Use np.load to load the npz file
data = np.load(npz_path)
```

---

## avg\_velocity.py

avg\_velocity.py contains:

- The `all_velocity` function, which calculates the average velocities between consecutive coordinates and stores them in a vector.
- The `plot_histogram` function, which creates and plots a histogram showing the frequency of each average velocity.

## Example Usage

### Usage

Use the `all_velocity` function to compute average velocities and the `plot_histogram` function to visualize the distribution of these velocities.

```
# Import necessary tools
import pytest
import numpy as np
import os
import pandas as pd
import matplotlib.pyplot as plt
from missing_data_dev.max_velocity.avg_velocity import calc_velocity,
avg_velocity, all_velocity, plot_histogram
```

```

# Create and declare the directory for all the CSV files
#direct_path =
"/Users/ibrahimrahat/Documents/GitHub/daphnia/data/table_data"
relative_path = "data/table_data"
all_files = os.listdir(relative_path)

# Create an empty dataframe
dataframes = []


# Loop through and have pandas read each CSV file
for file in all_files:
    file_path = os.path.join(relative_path, file)
    df = pd.read_csv(file_path)
    dataframes.append(df)

# Use the function all_velocity and store all the velocities into a
variable
all_velo = all_velocity(dataframes)

# Plot the velocities in a histogram
plot_histogram(all_velo)

```

```
In [ ]: display(Image.open('output_images/histogram of daphnia velocities.png'))
```





## Visualizing an Edge Discontinuity

- **Original Data:** Displays synthetic daphnia movement with edge discontinuities.
- **After Cleaning:** Data cleaned using TRexCleaner to address discontinuities.
- **After Imputation:** Data imputed with TRexImputer to smooth out gaps and improve continuity.

```
In [ ]: display(Image.open('output_images/side by side discontinuity example_highlighted.png'))
```

